



Sydney: 4 Carlotta Street, Artarmon, NSW 2064  
P: 02 8437 4600 F: 02 9437 4351

Melbourne: 1 / 35 Lakewood Blvd, Carrum Downs, VIC 3201  
P: 03 8773 4900 F: 03 9770 8241

## 24" Round Beret Headcover White Specification Sheet

**Item Code: 250311**

### Uses

A disposable Polypropylene headcover to protect the integrity of sensitive areas from contamination of scalp debris. (Hair, skin and dirt particles etc.)

Protective in all clinical areas requiring cross-contamination control.

### Features

Product Name:	Non-woven Circular Caps
Description:	Non-woven Circular, bouffant, beret caps
Application:	Dust protection
Colour:	White
Sizes:	24 inch
Fabric / Material:	Polypropylene Microfilament construct prevents material fibers and lint shedding. Soft and lightweight, breathable, air permeable. Latex Free
Fabric Weight:	14 gsm
Packaging Configuration:	1000 pcs per carton
Standards Conformity:	NA
MSDS:	NA

