

## 2189 GOWN

# SPECIFICATION

### FABRIC & CODE

ULTRASHIELD (US) ..... GDUS2189  
100% polyester with ECF stripe

### PROCESS | CLASS

Potentially Contaminated Garments Process Class

### STERILISATION

Yes

### FEATURES

- Knitted cuffs
- Clip closure at neck
- Tie at waist
- Suitable to wear when administering cytotoxic drugs

### SEAMS

Overlocked seams

### SIZES


SM - XL

### COLOUR

US Purple

### SUPPLY

Rental & Purchase



ALL SIZING SHOWN IN CENTIMETERS (CM) | ADDITIONAL SIZES AVAILABLE ON REQUEST

2189 GOWN		XSM	SM	M	L	XL
A	CHEST	-	127.0	147.3	157.5	172.7
B	SLEEVE	-	62.2	64.8	69.9	74.9
C	LENGTH	-	108.0	110.5	110.5	111.8