



6988 OVERSHOE

SPECIFICATION



FABRIC & CODE

ULTRASHIELD (US) GFUS6988

100% polyester with ECF stripe

DURACLEAN (PE)..... GFPE6988

100% polyester with ECF stripe

PROCESS | CLASS

US Potentially Contaminated Garments Process Class

PE Process Class B, D

STERILISATION

Yes

FEATURES

- Suitable to wear when dealing with cytotoxic drugs (UltraShield only)
- Aussie Tough sole

SEAMS

Overlocked seams

SIZES

XSM - 2XL

COLOUR

US Lemon, Green, Purple

PE Blue, Navy, Lemon, White

SUPPLY

Rental & Purchase

ALL SIZING SHOWN IN CENTIMETERS (CM) | ADDITIONAL SIZES AVAILABLE ON REQUEST

6988 OVERSHOE	XSM	SM	M	L	XL	2XL
A HEIGHT	47.0	47.0	50.8	50.8	57.2	57.2
B BOOT WIDTH DIAMETER	48.3	48.3	48.3	53.3	53.3	55.9
C LENGTH	26.0	28.0	30.0	32.0	34.0	36.6
D WIDTH	10.0	10.4	11.0	11.2	13.2	14.5