

# Klercide Neutral Detergent Unit Dose Concentrate

(3078810, 3078800)

## Technical Specification

PHYSICAL PROPERTIES	VALUE
Product Composition	Klercide Neutral Detergent Concentrate is a concentrated, sterile, pH neutral, non-ionic detergent with low residues.
Colour	Colourless
Odour	Slight
Clarity	Clear
Specific Gravity @ 20°C	0.990 – 1.020g/cm <sup>3</sup>
pH @ 20°C	5.5 – 7.0
Production	Klercide Neutral Detergent Concentrate is filtered to 0.2 micron and filled in a Grade C (ISO Class 7) cleanroom.
Packaging	All products are cleanroom packed, double bagged and supplied in 100ml bottles
Sterility	Klercide Neutral Detergent Concentrate is sterilised by gamma irradiation at not less than 25kGy.

TS052F-1