



T +44 (0) 2920 854 390
F +44 (0) 2920 854 391

BRUNEL WAY, BAGLAN ENERGY PARK,
NEATH, SA11 2GA, UK
www.ecolabcc.com
info@ecolab.com

Klercide Neutral Detergent

(3078770, 3078790, 3078780)

Technical Specification

PHYSICAL PROPERTIES	VALUE
Product Composition	Klercide Neutral Detergent is a sterile, pH neutral, non-ionic detergent with low residues.
Colour	Colourless
Odour	Faint
Clarity	Clear
Specific Gravity @ 20°C	0.990 – 1.010g/cm ³
pH @ 20°C	6.5 – 7.5
Production	Klercide Neutral Detergent is filtered to 0.2 micron and filled in a Grade C (ISO Class 7) cleanroom.
Packaging	All products are cleanroom packed and double bagged. Product code 3078770 is supplied in 5L capped bottles; product codes 3078790 and 3078780 are supplied in 1L trigger sprays. All trigger sprays are manufactured with the SteriShield Delivery System™
Sterility	Klercide Neutral Detergent is sterilised by gamma irradiation at not less than 25kGy.

TS049E-1