



T +44 (0) 2920 854 390
F +44 (0) 2920 854 391

BRUNEL WAY, BAGLAN ENERGY PARK,
NEATH, SA11 2GA, UK
www.ecolabcc.com
info@ecolab.com

Klercide 70/30 Denatured Ethanol (Aerosol) (3078380, 3078390)

Technical Specification

PHYSICAL PROPERTIES	VALUE
Product Composition	Klercide 70/30 Denatured Ethanol consists of a v/v formulation of 70% Ethanol blended with 30% Water for Injection.
Colour	Colourless
Odour	Alcohol-like
Clarity	Clear
Specific Gravity @ 20°C	0.850 – 0.890g/cm ³
Can Pressure	9-12 bar (0.2 micron filtered, oxygen-free, nitrogen propellant)
Endotoxin Limit	<0.25 EU/ml
Production	Klercide 70/30 Denatured Ethanol is filtered to 0.2 micron and filled into aluminium aerosol cans in a Grade C cleanroom.
Packaging	All products are cleanroom packed and triple bagged. This product is supplied in a 360ml aerosol can.
Sterility	Klercide 70/30 Denatured Ethanol is sterilised by gamma irradiation at not less than 25kGy.

TS034E-1