



T +44 (0) 2920 854 390
F +44 (0) 2920 854 391

BRUNEL WAY, BAGLAN ENERGY PARK,
NEATH, SA11 2GA, UK
www.ecolabcc.com
info@ecolab.com

Klercide 70/30 Denatured Ethanol

(3078260, 3078300, 3078290, 3078280)

Technical Specification

PHYSICAL PROPERTIES	VALUE
Product Composition	Klercide 70/30 Denatured Ethanol consists of a v/v formulation of 70% Denatured Ethanol blended with 30% deionised water
Colour	Colourless
Odour	Alcohol-like
Clarity	Clear
Specific Gravity @ 20°C	0.850 – 0.890g/cm ³
Refractive Index	19 – 21
Production	Klercide 70/30 Denatured Ethanol is filtered to 0.2 micron and bottled in a Grade C (ISO Class 7) cleanroom
Packaging	<p>All products are cleanroom packed and double bagged.</p> <p>Product code 3078280 is supplied in a 5L capped bottle; product code 3078290 is supplied as a 4L trigger spray; product code 3078300 is supplied as a 1L trigger spray; 3078260 is supplied as a 0.5L trigger spray.</p> <p>All trigger sprays are manufactured with the SteriShield Delivery System™</p>
Sterility	Klercide 70/30 Denatured Ethanol is sterilised by Gamma Irradiation at not less than 25kGy

TS021E-1