

Klercide 70/30 Denatured Ethanol

(3078340, 3078350, 3078360, 3078370)

Technical Specification

PHYSICAL PROPERTIES	VALUE
Product Composition	Klercide 70/30 Denatured Ethanol consists of a v/v formulation of 70% Denatured Ethanol blended with 30% Water for Injection
Colour	Colourless
Odour	Alcohol-like
Clarity	Clear
Specific Gravity @ 20°C	0.850 – 0.890g/cm ³
Refractive Index	19 – 21
Endotoxin Result	<0.25 EU/ml
Production	Klercide 70/30 Denatured Ethanol is filtered to 0.2 micron and bottled in a Grade C (ISO Class 7) cleanroom
Packaging	<p>All products are cleanroom packed and triple bagged.</p> <p>Product code 3078340 is supplied in a 5L capped bottle; product code 3078350 is supplied as a 4L trigger spray; product code 3078360 is supplied as a 1L trigger spray; 3078370 is supplied as a 0.5L trigger spray.</p> <p>All trigger sprays are manufactured with the SteriShield Delivery System™</p>
Sterility	Klercide 70/30 Denatured Ethanol is sterilised by Gamma Irradiation at not less than 25kGy

TS015E-1