

SINGLE LINEN SKIP

STAINLESS STEEL



FEATURES

- ▲ Trolley to suit Ironbach Bags
- ▲ Tubular stainless steel frame
- ▲ Overall dimensions: 460mmW x 540mmD x 870mmH
- ▲ Bag holder opening/ Internal dimensions: 395mmW x 355mmD
- ▲ 75mm Castors

ORDERING INFORMATION

CODE	DESCRIPTION	UNIT
DS1520SS	TUBULAR STAINLESS STEEL LINEN SKIP	EACH