



CSTMASK50PK

3-PLY SURGICAL NON-STERILE FACE MASK

SPECIFICATION SHEET

ITEM CODE

CSTMASK50PK

PRODUCT NAME:	3-PLY SURGICAL NON-STERILE FACE MASK
OUTER LAYER COMPONENT:	100% BLUE SBPP 23GSM
INNER LAYER COMPONENT:	100% WHITE SBPP 23GSM
FILTRATE COMPONENT:	PP WHITE FILTRATE PAPER (MELTBLOWN) 22GSM
ELASTIC:	NO LATEX ROUND WHITE ELASTIC
MASK LENGTH:	17.5CM
MASK WIDTH:	BEFORE PLEAT 17CM AFTER PLEAT 9.5CM
ALUMINIUM LENGTH:	10.5CM
RELAXED EARLOOP:	16CM
EARLOOP ELASTIC WIDTH:	0.6CM
COLOUR:	BLUE
PACKAGING CONFIGURATION:	50 PIECES PER PACK

FOR MORE INFORMATION PHONE 02 9851 4600 OR VISIT OUR WEBSITE www.cleanroom.com.au

TESTING

TEST	TEST STANDARD	STANDARD	TEST DATE	VALUE
Bacterial Filtration Efficiency (BFE)	EN14683:2019 (AnnexeB)	>98% (Type IIR)	JUNE, 18 2020	99.9%min
Breathability (DELTA P)	EN14683:2019 (AnnexeC) +AC: 2019	<60Pa/cm ² (TypeIIR)	JUNE, 18 2020	48max
Resistance	ISO22609:2004	≥16kPa	JUNE, 18 2020	Compliant
Intracutaneous irritation test	ISO10993-10	No irritation	JUNE, 18 2020	No irritation
Cytotoxicity	ISO10993-5	No cytotoxicity	JUNE, 18 2020	No cytotoxicity
Sensitisation	ISO10993-10 et-10A	Topical application: no sensitisation	JUNE, 18 2020	Topical application: no sensitisation
Sensitisation	ISO10993-10 et-10A	Intradermal injection: no sensitisation	JUNE, 18 2020	Intradermal injection: no sensitisation
Microbial cleanliness CFU/g	NF11737:2018	≤30 (Type IIR)	JUNE, 18 2020	13-18 CFU/g
Using time: BFE+Delta P	EN14683:2014 (AnnexeB&C)	>98%(Type IIR) <49Pa/cm ²	JUNE, 18 2020	8h:>99.9% 8h:<37,00Pa/cm ²