

Klercide™ Sterile 100% Polyester Isolator Cleaning Tool Pad



Klercide Sterile 100% Polyester Isolator Cleaning Tool Pads are designed to be used with the Klercide Isolator Cleaning Tool. They are for single use and supplied both gamma irradiated and double bagged. They are manufactured from 100% polyester.

PHYSICAL PROPERTIES	VALUE
Material Composition	<p><i>External Surface:</i> Double layer quilted 100% polyester, cleanroom laundered knitted fabric with laser sealed edges</p> <p><i>Internal Surface:</i> Single layer 100% polyester, cleanroom laundered knitted fabric with laser sealed edges</p>
Colour	White

PERFORMANCE DATA FOR EXTERNAL SURFACE OF THE CLEANING PAD		
Basic Weight	253g/m ² ±5%	TAPPI T-410
Absorbency		
Intrinsic	3.8ml/g	IES-RP-CC004.3
Extrinsic	967g/m ²	
Rate	1.0 sec	
Chemical Compatibility	Contains no chemical binders and is not affected by most common solvents & chemicals	

EXTRACTABLES	VALUE	PROCEDURE
Particles (>0.5µm)	0.92 x 10 ³ /cm ²	IES-RP.CC004.3
Fibres (>100µm)	0.14 fibers/cm ²	IES-RP.CC004.3
Non Volatile Residues		
Deionised Water	2.05 µg/cm ²	IES-RP.CC004.3
Isopropyl Alcohol	0.687 µg/cm ²	
Extractable Ions in DI Water		IES-RP.CC004.3
Chlorine (Cl ⁻)	None Detected	
Sulfate (SO ₄ ²⁻)	0.0282 µg/cm ²	
Nitrate (NO ₃ ⁻)	None Detected	
Sodium (Na ⁺)	0.0269 µg/cm ²	
Calcium (Ca ⁺)	0.0166 µg/cm ²	
Potassium (K ⁺)	None Detected	
Phosphate (PO ₄ ³⁻)	None Detected	

PACKAGING INFORMATION	
Packaging	Individually double bagged, 20 packs per case
Sterility	Sterile: Gamma irradiated at not less than 25kGy

PRODUCT ORDERING INFORMATION		
FORMAT	UNIT OF SALE	ORDER CODE
100% Polyester Isolator Cleaning Pads	20 pads	3038210

*Technical data represented in this table are typical values at the time of publication and within tolerances.

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

PLEASE SPEAK TO YOUR ECOLAB ACCOUNT MANAGER FOR FURTHER INFORMATION OR CALL
02920 854 390 (UK) +44 2920 854 395 (EXPORT)

WORLDWIDE HEADQUARTERS 1 Ecolab Place, St. Paul, MN 55102, USA www.ecolab.com/lifesciences
EUROPE HEADQUARTERS Richtistr. 7, 8304 Wallisellen, Switzerland www.ecolablifesciences.com

