



CRG MICRON™

**GAMMA IRRADIATED NON-WOVEN CLEANROOM WIPES
ISO CLASS 5-8 CLEANROOM COMPATIBLE**

Sterile Durx wipes (packs of 5). These changes are the result of 12 months of extensive validation, and enable CRG to offer cost reductions in combination with quality improvements



FEATURES

- High absorbency
- Low particles and extractables
- Validated sterile as per ISO 11137-2:2006
- Easy tear tabs
- Packaged in ISO Class 4 cleanroom
- Cleanroom compatible LDPE packaging
- IPA resistant printing
- Full traceability
- A Berkshire product manufactured in the UK and shipped to Ansell Malaysia to be re-packaged and sterilised (gamma) under the brand CRG Micron

INDUSTRY APPLICATIONS

- ✓ Pharmaceutical
- ✓ Medical Device
- ✓ Semi Conductor
- ✓ Electronics
- ✓ Pharmacy
- ✓ Food

ORDERING INFORMATION

OLD CODE	NEW CODE	SIZE	QUANTITY
0670S	CRGM0670	30 x 30 cm	5 WIPES PER INNER PK 150 WIPES PER OUTER PK
0770S	CRGM0770	30 x 30 cm	10 OUTER PKs PER CTN

FOR MORE INFORMATION PHONE 02 9851 4600 OR VISIT OUR WEBSITE www.cleanroom.com.au